

Viewing Patient History

On the Patient History tab, you can view a summary of the interactions the patient has experienced in the Pathways Portal tool. You can view details that help you understand where and when changes took place. For example, in the Alert History section, all of the alerts associated with the patient appear based on the status or severity. You can also export this list if needed for further investigation or record keeping.

Viewing Patient History

To view patient history:

1. On the **Patients** page, click in the patient header (profile area).

The Patient Profile page appears.

2. Click the **Audit** tab.


The Patient History page appears.


3. You can view the following types of patient history:

If you want to	Then
View all of the biometric data readings collected over a specific time period	Click Reading History .
View all of the patient's responses to their pathways questions over a specific time period	Click Pathway History .
View patient alerts based on status (accepted or rejected) or severity (high-red alerts, medium-yellow alerts).	Click Alert History .
View interactions the care team member completed in the patient's profile, such as adding a care plan or changing a minimum parameter for a biometric. This section can also contain actions the patient did or the system did automatically.	Click Profile History .

Expanding and Collapsing Sections

To expand or collapse history sections:


If you want to	Then
Open all of the patient history sections	Click  Expand All

If you want to	Then
Close all of patient history sections	Click 

Exporting a Patient History Report

To export a patient history record:

In the history section you want to export, click , and then click **MS Excel**.

 **Note:** You can export a Reading History report or a Pathway History report. Alert history for biometrics and pathway questions appear in these reports rather than in a separate report.

The report downloads to your computer in .xls format.

Biometrics Report													
Facility:		Vivify General											
Date Range:		12/30/2016-3/30/2017											
Patient Name:		Kent, Martha											
Encounter ID:		10											
Date/Time	Alert	Method	Device	BP	HR	SpO2	Wt	BS	Temp	PEF	FEV1	FEV 1/6	Alert Message
03/20/2017 12:42		Manual	Spiro									0.64	
03/20/2017 12:42		Manual	Spiro								2.64		
03/20/2017 12:42		Manual	Spiro							778			
03/20/2017 12:42		Manual	Temp						98.0				
03/20/2017 12:42	Medium	Manual	BP	127/102									BP 127/102 mm Hg (Exceeding or equal
03/20/2017 12:42		Manual	Wt				155.2						
03/20/2017 12:42		Manual	SpO2			96							
03/20/2017 12:42		Manual	SpO2		63								
03/20/2017 12:42		Manual	BS					95					
03/19/2017 12:42		Manual	Spiro									0.80	
03/19/2017 12:42		Manual	Spiro								4.68		
03/19/2017 12:42		Manual	Spiro							553			
03/19/2017 12:42	Medium	Manual	Temp						92.1				Temperature 92.1°F (Lower than or eq
03/19/2017 12:42		Manual	BP	141/76									
03/19/2017 12:42	Medium	Manual	Wt				155.4						Weight 155.4 lbs (3.4 lbs (2.14%) weigh
03/19/2017 12:42	Medium	Manual	SpO2			94							SpO2 94% (Lower than or equal to Aler
03/19/2017 12:42		Manual	SpO2		94								
03/19/2017 12:42	Medium	Manual	BS					127					Blood sugar 127 mg/dL (Exceeding or e